

# Information for powered wheelchair users

We realise that your powered wheelchair use might still be low, but it is important to remember to charge your batteries so your powered wheelchair is ready to use when needed.

Please follow these guidelines:



Powered wheelchairs should be charged every day, even after a small amount of use



Charge your wheelchair on a regular basis, even if you are not using it



Always make sure that when you turn your wheelchair on, the battery indicator shows full green lights



Refer to the manufacturer's instructions provided with your wheelchair for additional information. (This will have been supplied to you when your powered wheelchair was handed over)

**If you have any questions, please free to call us on:**

**York:  01904 654052**